Easy Telecart Guide





Quick Reference Guide Practitioner Cart

To Turn On Your Practitioner Cart

- To switch the Practitioner Cart on press the **Power**Button on the front of the Cart by the Battery
 indicator. You may need to also turn the HDX
 System on as well. To do this press the Power
 button on the HDX Remote control.
- The monitor may also need to be switched on.
- The system will take approximately one minute to fully load. The screen will then display the following screen:



Make a Call

Manually dialling the number:

- To enter the number you wish to dial manually, first select Place a Call from the VC Homepage.
- From the Place a Call screen, enter the number you wish to dial in the available field using your remote control.
- When you have entered the number and are ready to make your call, press the Call button (left of the directory button).

Dialling from the directory

- If the site you wish to dial to is in your directory press the **Directory** button in the centre of the remote control.
- Using the **Arrow** keys you can move through your directory list.
- When the specific number you wish to call is highlighted in Yellow press the Call button.
- The call will then try to connect.

Practitioner Cart Power



Flashing green lights:

The battery is charging.

All green Lights on:

The Battery is fully charged.

3 Green Lights:

30%-50% of Power left. An alarm will Beep to alert you.

Flashing Amber Light: 30% or less power in the cart. An alarm will continually

cart. An alarm will continually beep to warn you. Plug in the power cord in ASAP

No lights:

There is no power. You will have to plug in the power cord.



Aiming the HDX Camera

- To move and aim the camera both in and out of a call, press the Near button underneath the downward Arrow Key on the remote control. This will give you control of your camera.
- Using the Arrow keys you can pan the camera left / right and tilt it up / down.
- Using the Zoom buttons (to the right of the Arrow keys) you can zoom the camera in / out.

To Hang Up A Call

 Press the Hang-Up button, located to right of the Directory button.